

Product Name :
Filtration Assembly

Product Code :
NLAB-SCIENTIFICEQ18010



Description :

Filtration Assembly

Technical Specification :

Kit comprises of: All Glass Membrane Filter Holder, Borosilicate Glass, 47mm, with Sintered disc for membrane support and aluminum 'duck' clamp. Vacuum Pump, Compact and efficient.

Silicon Rubber Tubing 6x9mm -1 mtr.

Pack of 100 Membranes, 47mm.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003