

**Product Name :**  
Filtration Assembly

**Product Code :**  
NLAB-SCIENTIFICEQ18010



**Description :**

Filtration Assembly

**Technical Specification :**

**Kit comprises of:** All Glass Membrane Filter Holder, Borosilicate Glass, 47mm, with Sintered disc for membrane support and aluminum 'duck' clamp. Vacuum Pump, Compact and efficient.

Silicon Rubber Tubing 6x9mm -1 mtr.

Pack of 100 Membranes, 47mm.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003