

**Product Name :**  
Fisher Clamp

**Product Code :**  
NLAB-EDUSCHOOLEQ127001



**Description :**

Fisher Clamp

**Technical Specification :**

These Clamps are moulded in Polypropylene making it completely corrosion-proof.

The spring loaded jack in the middle provides a vice like grip to the burette and prevents it from slipping.

Two tapered wedges are provided to ensure firm grip on rods having dia 1/2" or 3/8".

Packing

Fisher Clamp - Single

61810r Clamp - Double

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003