

Product Name :
Fixed Voltage Regulator Trainer

Product Code :
NLAB-ELECTRICALAB210026



Description :

Fixed Voltage Regulator Trainer

Technical Specification :

SCOPE OF LEARNING:

Trainer used to study Voltage Regulation of

FIXED VOLATAGE Positive

FIXED VOLATAGE Negative

TECHNICAL SPECIFICATIONS:

Analog Meters

Voltmeter 30V

Ammeter 250mA

Power Supplies:

12-0-12 AC Isolated Power Supply.

Operated on Mains power 230V, 50Hz +10%

Components are mounted on the panels are:

Load Resistor using Band switch

Variable Resistor

SALIENT FEATURES:

Front panel built with high class insulated Printed Circuit Board sheet with well printed circuits and symbols.

Fuse for Short Circuit protection

Instruction manual.

Connections are brought out through 4mm Colored Sockets.

Patch Cords 4mm.

The trainer is housed in ABS Plastic cabinet.

Size of the trainer set 12"x8"

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003