

Product Name : Fixit Lever	Product Code : NLAB-LABORATORYEM77025
Description :	
Fixit Lever	
Technical Specification :	
<p>The levers can be fixed to the spindle in the center of which a vulcanite roller having 3 pin set at 90°on the central spindle.</p> <p>On one of these pins lever fits.</p> <p>Sliding through the frame is a rod carrying a spring that can be caused to bear against the roller its turning while other adjustments are made.</p> <p>Isotonic Lever:- with cam in place of after load screw without hook and weights.</p>	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003