

Product Name :
Flexometer Testing

Product Code :
NLAB-TESTINGLABM101001



Description :

Flexometer Testing

Technical Specification :

Smooth operation

Precisely engineered

High functionality

Angle of oscillation of grips $22.5^\circ \pm 0.5^\circ$

Flexing frequency 100 ± 5 cycles/minute

Dimensions of test specimens 70×45 mm

Maximum number of test specimen that can be tested at one time Six, twelve, or twenty-four

Motor $\frac{1}{4}$ HP single-phase 230 volts AC

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003