

Product Name :
Float Rack

Product Code :
NLAB-EDUSCHOOLEQ135001



Description :

Float Rack

Technical Specification :

Float rack molded in polypropylene is quite compact.

Autoclave able and can hold 16 micro centrifuge tubes for incubation in water baths.

Easy to load, there are essential in microbiology & immunology labs.

The surface of the rack has raised alpha-numeric index for easy identification and has a central rod handle for convenient lifting.

Material:- Polypropylene

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003