

**Product Name :**  
Float Rack

**Product Code :**  
NLAB-EDUSCHOOLEQ135001



**Description :**

Float Rack

**Technical Specification :**

Float rack molded in polypropylene is quite compact.

Autoclave able and can hold 16 micro centrifuge tubes for incubation in water baths.

Easy to load, there are essential in microbiology & immunology labs.

The surface of the rack has raised alpha-numeric index for easy identification and has a central rod handle for convenient lifting.

**Material:-** Polypropylene

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003