

**Product Name :**  
Flow Through Helical Coil India

**Product Code :**  
NLAB-ENGINEERINGLB36014



**Description :**

Flow Through Helical Coil

**Technical Specification :**

Helical Coil Material Stainless Steel Dia - 12"

No. of Turns 15

Coil Tube Dia OD 21mm, ID 17mm (Approx.)

Sump Tank Material Stainless Steel/PVC, Capacity 50 Ltrs

Water Circulation FHP Pump, Crompton/Sharp make

Manometer 2 Nos

Flow measurement Rotameter

Piping GI and PVC

Control panel comprising of Standard make on off switch, Mains Indicator etc

The whole set-up is well designed and arranged in a good quality painted structure

We are well-known manufacturers, OEM suppliers of Flow Through Helical Coil for Chemical Engineering Lab Equipments. Contact us for high quality Flow Through Helical Coil for Chemical Engineering Lab Equipments for schools lab, college lab, universities, research labs, various teaching and workshop training laboratories and industries in India.

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003