

Product Name :
 Fly Presses

Product Code :
 NLAB-MACHINELABWK26001



Description :

Fly Presses

Technical Specification :

Hand Fly Presses with two adjustable Vee guides provided with oil grooves for constant lubrication complete with 2 catch dogs with bolts and nuts. Fly wheel or handle with ball weights as desired (fly presses no. 0.1 & 2 are supplied with fly handle with one ball weight only or wheel as desired)

Standard Sizes

0.3/8"

3.3/2"

2"7/8"

3"1/2"

4"3/8"

5.3/8"

6.1/2"

7"5/8"

8.3/2"
9.3/4"
10.5/8"

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003