

Product Name :
Folding Machine

Product Code :
NLAB-MACHINELABWK24004



Description :

Folding Machine

Technical Specification :

Folding Machine .

CMT type folding machines are designed for easy and accurate clamping full length-wise bending and forming sheets to required angle correctly, angle bending stop has been fitted for correct repetition jobs.

REGISTRATION

48" **Width of sheet max.**

14 **Thickness of sheet max.**

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003