

**Product Name :**  
Food Technology Educational Lab Equipment

**Product Code :**  
NLAB-TECHNICALAB23001



**Description :**

Food Technology Educational Lab Equipment

**Technical Specification :**

Specifications:

Manual or interval timing control 0 - 10 minutes

Filter bag to retain stray particles

Power consumption: 3kw

Temperature range: 20 - 100°C ( $\pm 1^\circ\text{C}$ )

Maximum sample weight: 5kg

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

