

Product Name :
Food Technology Educational Lab Equipment

Product Code :
NLAB-TECHNICALAB23001



Description :

Food Technology Educational Lab Equipment

Technical Specification :

Specifications:

Manual or interval timing control 0 - 10 minutes

Filter bag to retain stray particles

Power consumption: 3kw

Temperature range: 20 - 100°C ($\pm 1^\circ\text{C}$)

Maximum sample weight: 5kg

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

