

Product Name : Forced Oscillation and Resonance	Product Code : NLAB-VOCATIONALTVET30010
Description :	
Forced Oscillation and Resonance	
Technical Specification :	
The system can be driven at various frequencies.	
Mass on the spring oscillates free or damped in water.	
Bright LED display of frequency.	
Easy to measure Amplitude.	
We are well-known manufacturers, OEM suppliers of Forced Oscillation and Resonance for Physics Workshop Labs. Contact us for high quality Forced Oscillation and Resonance for Physics Workshop Labs for schools lab, college lab, universities, research labs, various teaching and workshop training laboratories and industries in India.	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003