

Product Name :
Four Cylinder Balancing Machine

Product Code :
NLAB-TECHNICALAB35009

Description :

Four Cylinder Balancing Machine

Technical Specification :

Bench-top, to investigate vibration phenomena in a reciprocating multi-cylinder balancing machine.
Adjustable crankshaft and cylinder angles to study in-line, Vee and flat-four engines.
Adjustable piston mass and mass set.
Speed up to 1500rpm with all control gear and digital tachometer.
Frequency analyser to determine primary & secondary vibrations by accelerometer.
Interlocked safety canopy. 230V/50Hz.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003