

<b>Product Name :</b> Four Cylinder Balancing Machine	<b>Product Code :</b> NLAB-TECHNICALAB35009
<b>Description :</b>	
Four Cylinder Balancing Machine	
<b>Technical Specification :</b>	
<p>Bench-top, to investigate vibration phenomena in a reciprocating multi-cylinder balancing machine. Adjustable crankshaft and cylinder angles to study in-line, Vee and flat-four engines. Adjustable piston mass and mass set. Speed up to 1500rpm with all control gear and digital tachometer. Frequency analyser to determine primary &amp; secondary vibrations by accelerometer. Interlocked safety canopy. 230V/50Hz.</p>	

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003