

<b>Product Name :</b> Four Cylinder Four Stroke Petrol Engine with Rope Break Dynamometer Test Rig	<b>Product Code :</b> NLAB-TECHNICALAB36004
---	--



**Description :**

Four Cylinder Four Stroke Petrol Engine with Rope Break Dynamometer Test Rig

**Technical Specification :**

The setup consists of four cylinders, four stroke, Diesel engine connected to dynamometer for engine loading. The setup has stand-alone type independent panel box consisting of air box, fuel tank, manometer, fuel measuring unit. Engine jacket cooling water inlet, outlet and calorimeter temperature is displayed on temperature indicator. Flowmeters are provided for cooling water and calorimeter flow measurement.

**Experiment:**

Study of engine performance for power, efficiencies, fuel consumption, air fuel ratio and heat balance

**Utilities Required:**

Electric supply Provide 230 +/- 10 VAC, 50 Hz, single phase electric supply with proper earthing. (Neutral – Earth voltage less than 5 VAC)

5A, three pin socket with switch (2 No)

Water supply Continuous, clean and soft water supply @ 4000 LPH, at 10 m. head. Provide tap with 1" BSP

size connection

Space 3500Lx4000Wx2000H in mm

Drain Provide suitable drain arrangement (Drain pipe 65 NB/2.5" size)

Exhaust Provide suitable exhaust arrangement (Exhaust pipe 32 NB/1.25" size)

Fuel, oil Diesel @10 liter, Oil @ 3.5 lit. (15W40)

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003