

Product Name : Fried Ice Cream Machine	Product Code : NLAB-VOCATIONALTVET27063
Description :	
Fried Ice Cream Machine	
Technical Specification :	
Machine Body Material: Stainless Steel	
Rate Power Max: 1.5kW	
Voltage: 220 V	
Frequency: 50 Hz	
Dimension: 1100x600x800 mm	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003

