

Product Name :
FTIR Spectrometer

Product Code :
NLAB-SCIENTIFICEQ25007



Description :

FTIR Spectrometer

Technical Specification :

Spectral Range - 4000 to 400 cm⁻¹

Resolution 0.85cm⁻¹

Wavenumber

Scanning Speed for different applications

Signal to Noise Ratio (RMS value, at 2100cm⁻¹, resolution: 4cm⁻¹, detector: DTGS, 1 minute data collection)

Beam splitter

High efficiency Reflex

Detector

Data system computer

Software contains all routines needed for basic spectrometer operations, including library search, quantitation and spectrum export

IR Libraries include.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003