

Product Name :
Full Body Harness Edge 01

Product Code :
NLAB-SAFETYLABP20033



Description :

Full Body Harness Edge 01

Technical Specification :

Materials: Fray-proof dope-dyed polyester webbing Metal buckles of high strength alloy steel

Anchor Point: One dorsal D-ring, two chest attachment D-ring and two lateral D-ring.

Special Features: Ergonomic horizontal H-shaped thigh straps.

No discomfort or pressure on the groin.

Improved comfort in case of a fall.

Facilitate vertical and horizontal movements.

Breathable open-core padding for extra comfort throughout the neck, shoulder, back and thigh Easy no-tangle donning Built-in fall indicator for easy inspection

Stitching: High-strength polyester and computer-controlled stitching ensures consistency.

Adjustability: Size adjustment straps at chest, shoulders and thighs.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003