

Product Name :
Full Body Harness Tango I

Product Code :
NLAB-SAFETYLABP20035



Description :

Full Body Harness Tango I

Technical Specification :

Materials: Fray-proof dope-dyed polyester webbing Metal buckles of high strength alloy steel.

Anchor Point: One dorsal D-ring and two chest attachment textile loop.

Special Features: Ergonomic horizontal H-shaped thigh straps.

Facilitate vertical and horizontal movements.

No discomfort or pressure on the groin.

Improved comfort in case of a fall.

Built-in fall indicator for easy inspection Sub-pelvic straps stays low and comfortable.

Stitching: High-strength polyester and computer-controlled stitching ensures consistency.

Adjustability: Size adjustment straps at chest, shoulders and thighs.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003