

**Product Name :**  
Full Body Harness Tango I

**Product Code :**  
NLAB-SAFETYLABP20035



**Description :**

Full Body Harness Tango I

**Technical Specification :**

Materials: Fray-proof dope-dyed polyester webbing Metal buckles of high strength alloy steel.

Anchor Point: One dorsal D-ring and two chest attachment textile loop.

Special Features: Ergonomic horizontal H-shaped thigh straps.

Facilitate vertical and horizontal movements.

No discomfort or pressure on the groin.

Improved comfort in case of a fall.

Built-in fall indicator for easy inspection Sub-pelvic straps stays low and comfortable.

Stitching: High-strength polyester and computer-controlled stitching ensures consistency.

Adjustability: Size adjustment straps at chest, shoulders and thighs.

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003