

**Product Name :**  
Full Body Harnesses Di Electric Harness

**Product Code :**  
NLAB-SAFETYLABP20023



**Description :**

Full Body Harnesses Di Electric Harness

**Technical Specification :**

Materials: Fray proof dope : dyed polyester webbing. All metal components high class polymer coated for electrically insulated.

Anchor point: One dorsal D-ring and Two textile chest attachment loop and two lateral D-ring.

Special Features: Integrated work positioning waist pad ideally positioned sit-strap for extended comfort built-in dual color webbing for easy orientation.

Stitching: High-strength polyester stitching

Adjustability: Chest, Thigh & Shoulder

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

