

Product Name :
Full Body Harnesses Edge Lite Pro Harness

Product Code :
NLAB-SAFETYLABP20031



Description :

Full Body Harnesses Edge Lite Pro Harness

Technical Specification :

Anchor Point: One dorsal D-ring and one chest attachment D-ring.

Materials: Fray-proof dope-dyed polyester webbing Aluminum alloy quick buckles.

Special Features: Innovative lightweight design Extra comfort throughout the neck, shoulder and back
Breathable open-core padding technology Easy no-tangle donning Easy-to-wear with quick buckles Integrated
extension lanyard 500 mm long and D-ring Built-in fall indicator for easy inspection Sub-pelvic straps stays low
and comfortable.

Stitching: High-strength polyester and computer-controlled stitching ensures consistency.

Adjustability: Size adjustment straps at chest, shoulders and thighs.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

