

**Product Name :**  
Full Body Harnesses NC 1001 NM

**Product Code :**  
NLAB-SAFETYLABP20027



**Description :**

Full Body Harnesses NC 1001 NM

**Technical Specification :**

Stitching: High-strength polyester stitching.

Adjustability: Size adjustment straps only at chest and thighs.

Anchor Point: One dorsal D-ring and Two textile chest attachment loop.

Special Features: Ideally positioned sit-strap for extended comfort .

Built-in dual color webbing for easy orientation.

Materials: Fray-proof dope-dyed polyester webbing.

Metal parts covered with non-woven fabric pad.

Suitable for power sector work.

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003