

Product Name :
Full Body Harnesses NC 1103

Product Code :
NLAB-SAFETYLABP20046



Description :

Full Body Harnesses NC 1103

Technical Specification :

Anchor Point: One dorsal D-ring, Two Chest D-ring, Two Lateral D-ring.

Stitching: High-strength polyester stitching.

Adjustability: Size adjustment straps only at chest and thighs.

Special Features: Integrated work positioning waist pad Ideally positioned sit-strap for extended comfort Built-in dual color webbing for easy orientation.

Materials: Fray-proof dope-dyed polyester webbing Metal buckles of high strength alloy steel.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

