

Product Name :
Full Shoe Flex Tester Machines

Product Code :
NLAB-TESTINGLABM97001



Description :

Full Shoe Flex Tester Machines

Technical Specification :

Resistance of Footwear to Flexing

IS 8085 Method of Test for Footwear (Part 2) Footwear Performance Test, Stiffness Test for

Shanks; Tachometer Test for Cracking of Uppers; and Performance Test for Upper Fabrics, Coated Fabrics,

Sock Lining and other Lining Materials

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003