

Product Name :
Fully Automatic Autoclave

Product Code :
NLAB-LABORATORYEM47011



Description :

Fully Automatic Autoclave

Technical Specification :

Fully Automatic

Features:-

Automatic purging, Venting

Unique pedal free design

Flexible load sensor

Data logging of 2 sensor

Lower working height

Microprocessor based precise controlling

Low water level protection

Cooling can be either fast or slow

Self - diagnostic errors display

Fully automatic operation

User settable Multi programme cycle

We are leading manufacturers, suppliers of Fully Automatic Autoclave for Autoclave. Contact us to get high quality Fully Automatic Autoclave for Autoclave for schools, colleges, universities, research labs, laboratories and various industries.

Autoclave -

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003