

Product Name :
Fully Automatic Cryostat Microtome

Product Code :
NLAB-LABORATORYEM28012



Description :

Fully Automatic Cryostat Microtome

Technical Specification :

Chamber temperature as low as -40 deg C in less than 60 mins

Double Compressor with Peltier unit for fast freezing

Integrated Anti-roll device

UV-C light for chamber disinfectioning

Three speed sectioning

Auto/Manual Defrost

Programmable Start/Stop Time

Address: 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com