

Product Name :
Fully Automatic Double Column Horizontal Bandsaw

Product Code :
NLAB-MACHINELABWK18004



Description :

Fully Automatic Double Column Horizontal Bandsaw

Technical Specification :

STANDARD FEATURES :

Hydraulically controlled dual viceclamping system
Split front vice
Roller bearing / carbide blade guide system
Pre-set piece counter
Out of stock shut off
Automatic height control & rapid approach
Infinitely variable blade speeds by means of piv
Power driven band brush for chip removal
Hydraulic blade tension
Idler wheel motion detector with stalling & breakage of blade
Complete electricals including motors
Complete hydraulic system
Complete coolant system
Bimetal bandsaw blade
OPTIONAL ACCESSORIES :
Roller table 1500 mm
Multiple indexing attachment

Undle clamping attachment
Variable vice pressure

Working Capacity

Stroke 300mm
Stroke 120 to 1230
160mm

20 to 85 mm/in.

Motors :

Blade drive

Hydraulic pump

0.15HP pump

350 x 30 x 1.8mm.

Bar Feeding

Shuttle Bar Feed :

Single feed.

400 x 100mm

Machine weight

2150 x 2200 x 2160mm.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003