

**Product Name :**  
Fully Automatic Microtome

**Product Code :**  
NLAB-EDUSCHOOLEQ29005



**Description :**

Fully Automatic Microtome

**Technical Specification :**

Fully Automatic Microtome is a fully automatic microtome based on latest technology that functions on ball bearing guides and a microprocessor based control system with following specifications. It incorporates automatic memorization and recovery, a hand wheel which is lockable in any position, high volume trash tank, convenient specimen clamp, disposable embedding cassette clamp and paraffin clamp, safe alarm system, dormancy protection and advanced imported drive system.

Smartly located feather touch buttons can control all the functions effortlessly.

It can be used for semi-thin to thick sectioning of both soft and hard specimens.

**Technical Specifications:-**

Section Thickness Setting Range

Trimming Thickness Setting Range

Setting values

2-10 m increment step

10-20 m increment step

20-100 m increment step

|                       |
|-----------------------|
| Minimum setting value |
|-----------------------|

|                          |
|--------------------------|
| Vertical specimen stroke |
|--------------------------|

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003