

Product Name :
Fundamentals of Sedimentation India

Product Code :
NLAB-ENGINEERINGLB54021



Description :

Fundamentals of Sedimentation

Technical Specification :

- Sedimentation is often used to clarify suspensions.
- In the process, the solid particles move downwards in a liquid owing to their density.
- Using, the sedimentation processes in different suspensions can be investigated and compared.
- Five transparent cylindrical tanks are provided for the purpose.
- The suspensions are prepared in measuring cups, poured into the removable tanks, and mixed by shaking.
- The tanks are then mounted vertically on the experimental unit.
- To aid observation of the sedimentation process, the tanks are backlit.

We are well-known manufacturers, OEM suppliers of Fundamentals of Sedimentation for Applied Machines. Contact us for high quality Fundamentals of Sedimentation for Applied Machines for schools lab, college lab,

universities, research labs, various teaching and workshop training laboratories and industries in India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003