

Product Name :
Gas Cum Electric Autoclave

Product Code :
NLAB-LABORATORYEM47015



Description :

Gas Cum Electric Autoclave

Technical Specification :

Gas Cum Electric Autoclave -

Features:-

Biggest advantage where electricity is a major problem

Fully automatic operation

Energy efficient LPG operated Autoclave

Steam can be generated either by LPG or electricity

Specification:-

Capacity
300 x 500mm Dimension
Power Consumption
Operating/cm ² pressure
Operating 134°C Temperature

230V/50 HZ
Safety Valve
65kg Weight

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003