

**Product Name :**  
Gas Cum Electric Autoclave

**Product Code :**  
NLAB-LABORATORYEM47015



**Description :**

Gas Cum Electric Autoclave

**Technical Specification :**

Gas Cum Electric Autoclave -

**Features:-**

Biggest advantage where electricity is a major problem

Fully automatic operation

Energy efficient LPG operated Autoclave

Steam can be generated either by LPG or electricity

**Specification:-**

Capacity
300 x 500mm Dimension
Power Consumption
Operating/cm <sup>2</sup> pressure
Operating Temperature

230V/50 HZ
Safety Valve
65kg Weight

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003