

**Product Name :**  
Gauge Installation Case, Type MAK 1

**Product Code :**  
NLAB-EDUCATIONTVET15027



**Description :**

Gauge Installation Case, Type MAK 1

**Technical Specification :**

Thin and wide spatula

Scalpel blades

4H pencil

Cutting pliers

fillets

Accessories

2 rolls of transparent tape with dispensers

2 adhesive crepe paper with dispensers

1 roll of JG mylar tape (33 m)

1 teflon sheet (0.08x25x15 mm)

4 clamps N ° 1 (2.5 kg), opening 35 mm

2 clamps N ° 2 (2.5 kg), opening 60 mm

1 collet N ° 3 (7.5 kg), opening 80 mm

Metallic pressure plates

Glue

Glue M-200

Glue AE 10/15

Glue M-610

Soldering kit

Soldering iron with adjustable temperature, type Mark 5

Welding roll 361-A 20 R

Welding roll 450.20 S

Welding roll 570.28 R

AR flow kit

Assortment of CTF-AST relay lugs

Cable

30m coil (134 AWP)

30 m reel (126 DWV)

Reel of 30 m of cable, 3 flat conductors (326 DFV)

Reel of 15 m of cable, 1 teflon insulated conductor (130 FWT)

Protection products

2 bottles of M-coat A

1 bottle of M-coat B

1 bottle of M-coat C

1 bottle of M-coat D

1 bottle of M-coat G

1 tube of RTV 3145

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003