

Product Name : Gear Matic Hydraulic O.T. Table	Product Code : NLABMEDICAL-M2760002
Description : Gear Matic Hydraulic O.T. Table	
Technical Specification : The top segments are translucent for X-rays and enable taking pictures and mentioning a patient with a C-arm. The positions in five section table top are adjusted manually by mechanical drive. It's safe and robust mechanism prevents breakdown during an operation. Key Structure: With regular human strength, up-down table position can be fixed by stepping foot pedal downward The base is eccentric for trouble-free C-ARM movements. Leg & Head section are interchangeable Base & Column covered with stainless steel. Technical Constraints: Top Dimension: L 1905 X W 553 mm Height Adjustment: 787 mm -1042 mm Trandelenburg / Reverse: 30° / 25° Lateral Tilt: 20° / 20° Kidney Elevator: 150 mm Back Rest (up / Down): 80° / 25° Leg Rest (up / Down): 15° / 90°	

Head Rest (up / Down): 20° / 60°

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003