

Product Name :
Gear Matic Hydraulic O.T. Table

Product Code :
NLABMEDICAL-M2760002

Description :

Gear Matic Hydraulic O.T. Table

Technical Specification :

The top segments are translucent for X-rays and enable taking pictures and mentioning a patient with a C-arm.

The positions in five section table top are adjusted manually by mechanical drive.

It's safe and robust mechanism prevents breakdown during an operation.

Key Structure:

With regular human strength, up-down table position can be fixed by stepping foot pedal downward

The base is eccentric for trouble-free C-ARM movements.

Leg & Head section are interchangeable

Base & Column covered with stainless steel.

Technical Constraints:

Top Dimension: L 1905 X W 553 mm

Height Adjustment: 787 mm -1042 mm

Trandelenburg / Reverse: 30° / 25°

Lateral Tilt: 20° / 20°

Kidney Elevator: 150 mm

Back Rest (up / Down): 80° / 25°

Leg Rest (up / Down): 15° / 90°

Head Rest (up / Down): 20° / 60°

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003