

Product Name :
Gear Pump Cutaway

Product Code :
NLAB-AUTOMOTIVETVET21009



Description :

Gear Pump Cutaway

Technical Specification :

This cutaway model is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts and cross-sections. Many parts have been chromium-plated and galvanized for a longer life.

Automotive Technologies, Cutaways Models, Acrylic Static Cutaways Equipments for Automotive Training Lab

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003