


**Product Name :**  
Gearbox Apparatus

**Product Code :**  
NLAB-TECHNICALAB30016



**Description :**

Gearbox Apparatus

**Technical Specification :**

Measure of the angular displacement at the input and output of the gear box and comparison with the calculated values.

Calculation of the efficiency for different ratios according to the torque at entry.

**Technical specifications**

Fundamentally the gearbox consists of gear wheels of different sizes which may be engaged as required.

Simple three forward ratio and reverse sliding mesh box.

Pulleys fitted with protractors are attached to the input and output shafts so that the student may determine and verify velocity and torque ratios.

Experiment manual

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003