

**Product Name :**  
Geiger Muller Tube

**Product Code :**  
NLAB-EDUSCHOOLEQ146002



**Description :**

Geiger Muller Tube

**Technical Specification :**

A plastic cap is fitted over the vents for protection.

This is a halogen quenched metal tube with a thin end window, for detecting alpha, beta and gamma radiation and recommended as a general purpose tube for use with general school radioactive sources

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003