

Product Name :
Glass Microfiber Filter Paper Whatman

Product Code :
NLAB-SCIENTIFICIND13004



Description :

Glass Microfiber Filter Paper (Whatman)

Technical Specification :

45µm pore size or Cellulose Nitrate, suitable for water quality ;

membrane diameter 47mm pore size;

white in color

These depth filters combine fast flow rates with high loading capacity and the retention of very fine particle, extending into the sub-micron range.

Ideal for use in applications involving air filtration and for gravimetric analysis of volatile materials where ignition is involved

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

