

**Product Name :**  
Glass Microfiber Filter Paper Whatman

**Product Code :**  
NLAB-SCIENTIFICIND13004



**Description :**

Glass Microfiber Filter Paper ( Whatman)

**Technical Specification :**

45µm pore size or Cellulose Nitrate, suitable for water quality ;

membrane diameter 47mm pore size;

white in color

These depth filters combine fast flow rates with high loading capacity and the retention of very fine particle, extending into the sub-micron range.

Ideal for use in applications involving air filtration and for gravimetric analysis of volatile materials where ignition is involved

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

