

Product Name :
Glass Syringe

Product Code :
NLAB-EDUSCHOOLEQ46006



Description :

Glass Syringe

Technical Specification :

We carry the most comprehensible line of glass syringes with multi fit interchangeability and are supplied with non-sterilization condition.

The quality of our glass syringes is unbeatable quality.

These are available in various capacities 1ml / 2ml / 5ml / 10ml / 20ml / 30ml / 50ml / 100ml.

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments