

Product Name : Glass Tube For Physics Lab	Product Code : NLAB-EDUSCHOOLEQ1730194
Description : Glass Tube For Physics Lab	
Technical Specification : Used in laboratory, Heat resistant, Chemical resistant, Used for stirring and mixing. The Glass Tube for Physics Lab Both ends closed and rounded. Material: Borosilicate glass Dimensions : 300 x 15 mm.	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003