

Product Name :
Glue Spreader

Product Code :
NLAB-MACHINELABWK33030



Description :

Glue Spreader

Technical Specification :

Main rubber roll: Ø230 x 1400mm

Doctor roller: Ø188 x 1400mm

Lower roller: Ø230 x 1400mm

Clearance between rubber and rollers: 35mm

Feeding speed: 32 m/min

Must come complete with an air compressor with 3.7kW motor, which require 7 kg/cm air pressure of 1200 l/min free air.

Wood Workshop Portable Tools Glue Spreader Manufacturer, Supplier and Glue Spreader Exporter in India

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003