

Product Name :
Graduated Measuring Jars

Product Code :
NLAB-LBEMALAB230003



Description :

Graduated Measuring Jars

Technical Specification :

These are commonly used at site for measurement of Silt & Sand, to measure the quantities of liquids these are available in Ordinary Glass, Polypropylene Plastic, Borosilicate Glass, Soda Glass. Capacity available form 5ml to 2000ml.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003