


**Product Name :**  
Graduated Measuring Jars

**Product Code :**  
NLAB-LBEMALAB230003



**Description :**

Graduated Measuring Jars

**Technical Specification :**

These are commonly used at site for measurement of Silt & Sand, to measure the quantities of liquids these are available in Ordinary Glass, Polypropylene Plastic, Borosilicate Glass, Soda Glass. Capacity available form 5ml to 2000ml.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003