

Product Name :
Grains Moisture Meter

Product Code :
NLAB-COMPSLAB280004



Description :

Grains Moisture Meter

Technical Specification :

Accuracy: $\pm 0.2\%$

Material: Stainless Steel

Size: 31 X 18 X 19 Cm approx

Range: 8% to 40%

Power: 230 V

Display: Digital 3 Digits LED

Food Grains

Pulses

Oil Seeds

Vegetable Seeds,

Rice

Paddy

Wheat

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003