

**Product Name :**  
Grinding Machine and Accessories

**Product Code :**  
NLAB-MACHINELABWK34005



**Description :**

Grinding Machine and Accessories

**Technical Specification :**

Type: Double Head type pedestal grinder

Grinding Wheel: Wheel Size: Approx. Dia. 300mm

Wheel Grit: Approx. No.36 and No. 60

(one fine and one coarse)

Motor output: 2HP

Spindle Speed: Approx. 1500 rpm

Pedestal Stand

Coolant System: - coolant pump; coolant tray; coolant tank and fittings (must be fitted in the stand)

Tool Tray:- Must be equipped at each side

Spark Protector: Must be equipped at each side

Wheel Guard: Must be equipped at each side

Power Source: 220V 50Hz, Single phase

Accessories:

Standard Accessory set

Spare Grinding wheels

Grit No. 36 and No. 60, 3 pcs each

Operations manual

Standard Welding Equipment Laboratory Set Grinding Machine and Accessories  
Manufacturer, Supplier and Grinding Machine and Accessories Exporter in  
India

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003