

Product Name :
GTM Screen Stability Test

Product Code :
NLAB-CONCTLAB240003



Description :

GTM Screen Stability Test

Technical Specification :

10 Litre Bucket with lid

5 mm sieve 350 mm diameter

Sieve pan

Balance, 20 g minimum capacity 20 kg

Stopwatch

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003